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ARIA® e Eclipse™

Thais Flores

Clinical Applications Specialist LATAM

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Improving radiotherapy safety and efficiency with the customized ARIA oncology information system

Zhongfei Wang¹, Qinghui Yun², Changhao Liu¹, Xiaohuan Sun¹, Wei Wang¹, Yutian Yin¹, Feng Xiao¹, Lina Zhao¹

Affiliations + expand
PMID: 34421003 DOI:

Received: 14 February 2019 | Revised: 14 February 2019 | Accepted: 15 February 2019
DOI: 10.1002/acm2.12555

PARALLEL OPPOSED EDITORIAL

WILEY

Clinical practice workflow in Radiation Oncology should be highly standardized

INTERNATIONAL JOURNAL OF
RADIATION ONCOLOGY • BIOLOGY • PHYSICS ASTRO

POSTER VIEWING | VOLUME 99, ISSUE 2, SUPPLEMENT, E507, OCTOBER 01, 2017

The Importance of Standardization in Workflow Enhancement and Patient Safety

C.M. Leyrer • J. Kittel • T. Meier • ... J.H. Suh • R.D. Tendulkar • S.

DOI: <https://doi.org/10.1016/j.ijrobp.2017.06.1939> Check for updates



Received: 14 January 2021 | Revised: 21 May 2021 | Accepted: 28 June 2021
DOI: 10.1002/mp.15116

AAPM SCIENTIFIC REPORT

MEDICAL PHYSICS

Electronic charting of radiation therapy planning and treatment: Report of Task Group 262

James G. Mechalakos¹ | Sonja Dieterich² | Luis E. Fong de los Santos³ |
Sandra C. Fontenla¹ | Joseph Hanley⁴ | Vijay A. Harwalkar⁵ | Linda X. Hong⁶ |
Y. Jessica Huang⁷ | Grace Gwe-Ya Kim⁸ | Susan L. Richardson⁹ |
Steven G. Sutcliffe¹⁰ | Sridhar Yaddanapudi¹¹ | Sue Merkel¹² | Mark Parry¹³

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

Appointment Scheduling Care Path Registration Diagnosis Summary Encounters Chart QA Diário

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Hospital

Hospital principal Departamento principal
Las Vegas General Radiation Therapy

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Alertas

ARIA®

The screenshot displays the Varian ARIA software interface. A modal window titled "Alerta de paciente" is open, showing patient information for Jennifer Gilliano (US-STT 4656468). The "Alerta de" field is highlighted with a red box and contains the text: "Encaminhar paciente para revisão médica e avaliação de toxicidade". The modal also includes a "Data de" field set to 7/27/2022 and buttons for "OK" and "Cancelar".

The background interface shows a patient summary for Jennifer Gilliano, including contact information, hospital details, and a treatment plan. The treatment plan table is as follows:

ID do plano	ID Rx	Fração	Prescribe Details	Prescribe To	Dose Aplicada (cGy)	Dose Restante (cGy)	Primeiro Tratame...	Último Tratamento
C2 Pelvis	Pelvis 4FLD	Pelvis	3/25	Isocenter	540	3960	11/10/2015	11/11/2015
C2 Pelvis	Pelvis 4FLD	Pelvis	25	Isocenter				
C1	(nenhum)	Mama	15	Volume				
C1	4 Field		7/25		1260	3240	2/9/2009	2/17/2009
C1	2 Isocenters		0/10		0	2000		
C1	Delta couch		0/25		0	4500		

Below the treatment plan, there is a section for "Pontos de Referência" and a table for "Alertas do Paciente". The "Alertas do Paciente" table shows a new alert created on 7/27/2022 at 8:05:43 AM, with the text "Encaminhar paciente para revisão médica e avaliação de toxicidade". A blue arrow points to the "Criado por" field, which contains the name "cwalker".

Alertas

ARIA®

The screenshot displays the Varian ARIA console interface. At the top, it shows 'Medical Systems - 4D Console version 13.0', 'Clinac IX #3649 VAR_JEC Scale', and the date '27-Jul-2022'. The main window is titled 'varian Queue' and shows a patient named 'Gilliano Jennifer' with a photo and personal details: Idade: 48, Data de Nascimento: 21/02/1974, ID1: US-STT 4656468, ID2: 514896, Radioterapeuta: dvance, and Diagnóstico: ICD-10:D06.9 Carcinoma in situ: Cervix, unspecified.

A table titled 'PROGRAMAÇÃO DE HOJE PARA 21EX 1900' shows the patient's schedule. A blue arrow points to a warning icon in the 'Atividade' column for the entry at 10:00.

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Below the table, there are buttons for 'Remover Paciente de Hoje' and 'Atualizar'. To the right, there is a 'Curso' section with a list of courses and fields, and a 'Tratar' button. At the bottom right, there is a 'Plan' section with a table for 'Couch Vrt(cm)', 'Couch Lng(cm)', 'Couch Lat(cm)', 'Couch Rtn(deg)', and 'SSD(cm)', and a 'Clinac Console' section with buttons for 'RV Mode Up', 'Goto', 'Beam On', 'Beam Off', 'Delta Mode', 'Auto-Mode Up', and 'Always On Top'.

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Alertas

ARIA®

The screenshot displays the Varian ARIA console interface. At the top, it shows 'Medical Systems - 4D Console version 13.0', 'Clinac_TX #3649 VAR_JEC Scale', and the date '27-Jul-2022'. The main window is titled 'varian Queue' and shows a patient's profile for 'Gilliano Jennifer'. The patient's details include: Idade: 48, Data de nascimento: 21/02/1974, ID1: US-STT 4656468, ID2: 514896, Radioterapeuta: dvance, Diagnóstico: ICD-10:D06.9, and Carcinoma in situ: Cervix, unspecified. A table titled 'Alertas do Paciente' shows a single alert for '27 de jul de 2022' with the message 'Encaminhar paciente para revisão médica e avaliação de toxicidade'. A blue arrow points to the 'Aplicar' button at the bottom of the alert window. The right side of the interface shows treatment parameters: Couch Vrt(cm) 0.0, Couch Lng(cm) 0.0, Couch Lat(cm), Couch Rtn(deg), and SSD(cm). A 'Clinac Console' panel on the right contains buttons for 'RV Mode Up', 'Goto', 'Beam On', 'Beam Off', 'Delta Mode', 'Auto-Made Up', and 'Always On Top'. The bottom of the interface shows a 'Comp Mount' section with a password field and a 'Bols' section.

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Show Setup Note...

Tratar CQ UTM

RV Mode Up
 Goto
 Beam On
 Beam Off
 Delta Mode
☒ Auto-Made Up
☒ Always On Top

ARIA®

N MEDICAL SYSTEMS, INC. For immediate recipient's use only.

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
5.8 | Prescription

Users should take advantage of the capabilities for prescribing that are provided by the RO-EMR. ASTRO (American Society for Radiation Oncology) has provided guidance on items to include to improve standardization of dose prescriptions.³¹

Dose-volume constraints can be considered as part of the prescription or as a separate document of intent to ensure that treatment planning obtains the complete information to begin the plan, thus limiting the need for unnecessary communication or revision downstream. These constraints can be explicitly stated or standard department constraints can be referenced with explicit exceptions listed for the particular case. Any modification to constraints will result in a prescription modification if the document is approved. Institutional standards for personalized dose constraints alleviate this inefficiency by requiring explicit documentation of constraints which differ from the institutional norm.

Common failure modes associated with the prescription fall into three general categories:

1. Incomplete information and typographical errors
2. Mismatches between the prescription and the treatment plan
3. Changes to the prescription that are not communicated to the clinic

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


4.1.4.3. Minimizing Time Pressures

In order to avoid safety problems or quality lapses caused by rushing to meet unrealistic scheduling expectations, each practice should determine the appropriate time allocated for each step in the process of care. While the time needed to perform processes will vary between practices, **Table 4.1** provides a template, listing basic steps in the process. The goal of this effort is to avoid safety issues caused by time pressures, while satisfying the responsibility of the clinical team to set a course of action that will assure a timely, yet safe and accurate transition from patient clinical evaluation to treatment.

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Varian Medical Affairs

MyVarian

The image displays two overlapping screenshots of the Varian Medical Affairs website. The top screenshot shows the 'myvarian' portal with a navigation bar containing 'Documentation', 'Support', 'Training', 'Events', and 'Resources'. A dropdown menu for 'Varian Medical Affairs' lists 'Marketing Your Center', 'Service & Education', and 'Varian Medical Affairs'. The main content area features a blue banner with the text 'Do you have an idea to improve a Varian product?' and a call to action to visit the Product Ideas workspace. The bottom screenshot shows the main Varian website with a navigation bar including 'MyVarian', 'Medical Affairs', 'Investors', 'For Patients', and 'Careers'. The 'Medical Affairs' section is highlighted, and a 'Tools' dropdown menu is open, listing 'RapidPlan Model Directory', 'Protocol Clinical Goals', 'Example Dosimetric Scorecards', and 'Community Software'. The main content area on the bottom screenshot is titled 'VARIAN MEDICAL AFFAIRS MISSION:' and describes the goal of engaging with medical professionals to exchange clinical information and tools to beat cancer. The Varian logo and tagline 'A Siemens Healthineers Company' are visible in the bottom right corner of both screenshots.

English | English ▾ Varian.com thais.pires.flores@hotmail.com Contact Us My Account Log Out

myvarian Documentation Support Training Events Resources

Do you have an idea to improve a Varian product?
Visit the Product Ideas workspace to submit your idea and to vote for ideas already submitted.

Marketing Your Center
Service & Education
Varian Medical Affairs

varian
A Siemens Healthineers Company
A leap in cancer care. act.

MyVarian Medical Affairs Investors For Patients Careers

varian
A Siemens Healthineers Company

Why Varian Products Resources & Support About Varian Search

Medical Affairs Studies Literature Meetings Tools Contact

RapidPlan Model Directory
Protocol Clinical Goals
Example Dosimetric Scorecards
Community Software

Feedback

VARIAN MEDICAL AFFAIRS MISSION:

Engage with medical professionals to exchange meaningful clinical information, tools, scientific research, and expert opinion to beat cancer and other conditions utilizing radiation therapy. We strive to maximize the clinical potential of our products, research and development efforts and collaborative research partnerships with clinical investigators.

Through our relationships with clinicians and researchers:

varian
A Siemens Healthineers Company

Varian Medical Affairs

MyVarian

Publicly shared RapidPlan models

Please feel free to share your model via [PeerAnswers](#) and we will add it the next time this list is updated

Tx Site	Institution	PDF	Publication	Rx (max/nom)	Fraction	Protocol link	Share link	Modality	Version
Brain-HSWBv2	Varian MedAffairs	yes	Liu H et al.	30Gy	3Gy	NRG-CC001	VarianMedicalAffairs	Photon	15.6 + 17
Brain-HSWBv1	Varian MedAffairs	yes	Magliari V et al.	30Gy	3Gy	NRG-CC001	VarianMedicalAffairs	Photon	13.5
Prostate	University of Miami	no		80Gy (SIB)	2Gy	similar Rx referenced	Peer Answers	Photon	13.6.23
R Lung SBRT	UC San Diego	yes		56Gy	8Gy	similar Rx referenced	Peer Answers	Photon	13.6.23
Head-and-Neck	UC San Diego	yes		70Gy	2Gy	RTOG0615	Peer Answers	Photon	13.6.23
Prostatic Fossa	UC San Diego	yes		70.2Gy (SIB)	1.8Gy	RTOGP0126	Peer Answers	Photon	13.6.23
L Lung SBRT	UC San Diego	yes		56Gy	8Gy	similar Rx referenced	Peer Answers	Photon	13.6.23
Prostate SV	UC San Diego	yes		81Gy	1.8Gy	similar Rx referenced	VarianMarketplace	Photon	13.6.23
Prostate SV	University of Miami	no		80Gy (SIB)	2Gy	similar Rx referenced	VarianMarketplace	Photon	13.6.23
Prostate SV LN	University of Miami	no		80Gy (SIB)	2Gy	similar Rx referenced	VarianMarketplace	Photon	13.6.23
Esophagus	Varian MedAffairs	yes	Magliari A	50.4Gy	1.8Gy	RTOG1010	VarianMedicalAffairs	Photon	15.6
Esophagus	Varian MedAffairs	yes	Celic E et al.	41.4Gy	1.8Gy	see publication	VarianMedicalAffairs	Proton	16.1
Liver	Varian MedAffairs	yes	Cozzi L et al.	50Gy	10Gy	NRG-G1003	VarianMedicalAffairs	Proton	16.1
Head-and-Neck	Varian MedAffairs	no	Delaney A et al.	70Gy (SIB)	2.12Gy	RTOG0615	VarianMedicalAffairs	Proton	16.1
Head-and-Neck	Northampton General	no		65Gy (SIB)	2.2Gy	RTOG0234	Peer Answers	Photon	15.6.04
Cervical cancer	Royal Surrey County	yes	Hussein M et al.	50.4Gy	1.8Gy	INTERTECC, INTERLACE	Peer Answers	Photon	13.5.15
Head-and-Neck	Univ College London	no		66Gy (SIB)	2Gy	RTOG0234	Peer Answers	Photon	15.6.04
Anorectal	UC San Diego/OrbitRT	yes		62Gy (SIB)	2.21Gy	APHRODITE	ORBIT-RT	Photon	15.6
Breast APBI	UC San Diego/OrbitRT	no		30Gy	3Gy BID	NSABP B-39 RTOG0413	ORBIT-RT	Photon	16.1
Glioblastoma	UC San Diego/OrbitRT	no		60Gy (SIB)	2Gy	RTOG0825	ORBIT-RT	Photon	15.6
GYN multi-target	UC San Diego/OrbitRT	no		50.4Gy (SIB)	1.8Gy	RTOG0418	ORBIT-RT	Photon	15.6
GYN single target	UC San Diego/OrbitRT	no		50.4Gy	1.8Gy	RTOG0418	ORBIT-RT	Photon	15.6
Head-and-Neck	UC San Diego/OrbitRT	no		70Gy	2Gy	RTOG0615	ORBIT-RT	Photon	16.1

[Tools](#)

RapidPlan
Model
Directory

Protocol
Clinical
Goals

Example
Dosimetric
Scorecards

Community
Software

Varian Medical Affairs

MyVarian

Eclipse V16+ Clinical Goal Templates

Patient Site	Rx	Fraction	Protocol	Download Clinical Goal XML
GYN	50.4Gy	1.8Gy	NRG GY006	NRG-GY006_50.4Gy
Prostate	79.2Gy	1.8Gy	RTOG0126	RTOG0126_79.2Gy
Anus+LN	54Gy	1.8Gy	RTOG0529	RTOG0529_54Gy
Brain	60Gy	2Gy	NRG BN001	NRG-BN001_60Gy
Brain	40Gy	2.66Gy	NCIC CE.6	NCI-CCE6_40Gy
Esophagus	50Gy	2Gy	EA2174	EA2174_50Gy
Prostate+LN	70Gy	2.5Gy	NRG GU005	NRG-GU005_70Gy
Head-and-Neck	70Gy	2.12Gy	RTOG0615	RTOG0615_70Gy
Head-and-Neck	70Gy	2Gy	NRG HN004	NRG-HN004_70Gy
Rt Breast APBI	30Gy	6Gy	Florence trial	Florence_30Gy
Lt Breast APBI	30Gy	6Gy	Florence trial	Florence_30Gy
Rt Breast APBI	30Gy	6Gy	ACCEL Trial	ACCEL_30Gy
Lt Breast APBI	30Gy	6Gy	ACCEL Trial	ACCEL_30Gy
Rt Breast	48Gy	3.2Gy	RTOG1005	RTOG1005_48Gy
Lt Breast	48Gy	3.2Gy	RTOG1005	RTOG1005_48Gy
Lung	60Gy	2Gy	NRG LU004	NRG-LU004_60Gy
Lung	60Gy	4Gy	NRG LU004	NRG-LU004_60Gy
Pancreas	50.4Gy	1.8Gy	RTOG0848	RTOG0448_50.4Gy
Prostate Fossa	70.2Gy	1.8Gy	NRG GU006	NRG-GU006_70.2Gy
R Breast/CW+LN	50Gy	2Gy	NSABP51/RTOG1304	NSABP51-RTOG1304_50Gy
L Breast/CW+LN	50Gy	2Gy	NSABP51/RTOG1304	NSABP51-RTOG1304_50Gy
Rectum	50Gy	2Gy	UK IMRT Guidance	UK-IMRT-Guidance_50Gy
Rectum	25Gy	5Gy	UK IMRT Guidance	UK-IMRT-Guidance_25Gy
Liver SBRT	50Gy	10Gy	RTOG1112	RTOG1112_50Gy
Lung SBRT	34Gy	34Gy	RTOG0915	RTOG0915_34Gy
Lung SBRT	48Gy	12Gy	RTOG0915	RTOG0915_48Gy
Prostate SBRT	36.25Gy	7.25Gy	NRG GU005/RTOG0938	NRG-GU005/RTOG0938_36.25Gy

[Tools](#)

RapidPlan
Model
Directory

Protocol
Clinical
Goals

Example
Dosimetric
Scorecards

Community
Software

Varian Medical Affairs

MyVarian

Varian Innovation Center GitHub projects

These community co-developed and supported tools are hosted on the VIC GitHub instance under the [Varian Limited Use Software License Agreement](#). The below list of publicly available software projects / source code are not in any way validated or verified for clinical use by Varian. Always perform due diligence, testing and validation of any tools used in your clinic.

Project	Publication	Description
PlanScoreCard		Medical Affairs Applied Solutions ESAPI tool to create ScoreCards and score plans; in-metric Boolean/expans
DoseDynamicArcs (coming soon)		Medical Affairs Applied Solutions ESAPI tool to edit dose rate and gantry speed using various methods; enabl

Individual public GitHub clinical projects

These community developed and supported tools are hosted on the public GitHub source code repository, usually under open source MIT License agreement. The below lists of publically available software projects / source code are not in any way endorsed, validated, verified or otherwise tested by Varian. Always perform due diligence, testing and validation of any tools used in your clinic.

Project	Publication	Description
VMAT-TBI	Simiele E	ESAPI code used for VMAT TBI autoplanning
eFLASHBeamModeltoTPS	Rahman M et al.	Electron FLASH beam model
PlanCheck		A Varian Eclipse Scripting API plugin that checks DVH metrics, potential collisions, and plan/structure/dose p
AutoRing_SRS		Generates UAB rings (inner, middle, outer) from a PTV.
DVHAsync		Compares DVH metrics between two plans or plansums.

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

Received: 3 June 2021 | Revised: 14 September 2021 | Accepted: 30 October 2021
DOI: 10.1002/acm2.13479

RADIATION ONCOLOGY PHYSICS

JOURNAL OF APPLIED CLINICAL
MEDICAL PHYSICS

Varian ethos online adaptive radiotherapy for prostate cancer: Early results of contouring accuracy, treatment plan quality, and treatment time

Mikel Byrne^{1,2} | Ben Archibald-Heeren¹ | Yunfei Hu³ | Amy Teh¹ |
Rhea Beserminji¹ | Emma Cai¹ | Guilin Liu¹ | Angela Yates¹ | James Rijken⁴ |
Nick Collett¹ | Trent Aland⁵

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The moment of
Ethos™ therapy
is now.

Yesterday's belief in how cancer should be fought has changed forever. From this day forward, we're committed to helping you put the power of **Adaptive Intelligence™** to work in your fight against cancer so you can go beyond...

Beyond... complex adaptive treatments to a new, AI-driven therapy solution.

Beyond... what you thought possible today and toward a new belief—a new ethos—in what you'll be able to achieve, impact, and inspire on a global scale.

The moment of Ethos™ therapy is now.

Visit varian.com/ethos to learn more and experience Ethos™ therapy, an Adaptive Intelligence solution from Varian.

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